	<p>Innovative Optics LLC 6812 Hemlock LN Maple Grove MN 55369 Ph: 763-425-7789 Fax: 763-425-6689 Email: sales@innovativeoptics.com www.innovativeoptics.com</p>
---	---

**EU-Declaration of
 Conformity For: PT-CEIII
 Overmold Patient Goggle**

Manufacturer: Innovative Optics, LLC
 6812 Hemlock Ln Maple Grove MN 55369


Innovative Optics declares the part listed below is in conformance with the provisions of the PPE Regulation (EU) 2016/425 Directive 2001/95/EC and where such is the case, with the harmonized standard EN 207:2017

Model Number: PT-CEIII Overmold Patient Goggle

Description:

Stainless steel patient goggle.

- Our new technology allows clinicians to adjust interpupillary distance for different patients, while it is on the patient.
- The high-arching nose bridge pivots up and down to allow for tight procedures around the eyes. ackaging consists of the PT-CEIII Overmold Patient Goggle, one white silicone headstrap, and carrying case.
- This goggle can be sterilized in solution or autoclaved.

<p>> 315-1400 D LB6 + IR LB8 + M LB8Y INVO CE > 1400-11000 D LB4 + I LB5 INVO CE UKCA</p>	
--	---

The product is identical to the PPE which is subject of the EU-type-examination certificate no C6221INVO/R1

<p>DIN CERTCO Gesellschaft für Konformitätsbewertung mbH Alboinstr. 56 12103 BERLIN GERMANY Identification number of the Notified Body: 0196</p>	
---	--

The declaration of conformity is issued under the sole responsibility of the manufacturer

Signature: 

Thomas Barrows
 President

Date Issued: 7.30.21