



Innovative Optics  
 6812 Hemlock LN Maple Grove MN 55369  
 Ph: 763-425-7789 Fax: 763-425-6689  
 Email: [sales@innovativeoptics.com](mailto:sales@innovativeoptics.com)  
[www.innovativeoptics.com](http://www.innovativeoptics.com)

**Declaration of Conformity**  
**For: Renegade**

Manufacturer: Innovative Optics, Inc  
 6812 Hemlock Ln Maple Grove MN 55369

Innovative Optics declares the part listed below is in conformance with the provisions of the PPE Regulation 2016/425/EU of the European Parliament and of the council, Directive 2001/95/EC and where such is the case, with the harmonized standard EN 207:2017

**Model Number: Renegade**

Description:

The Renegade is a Laser Safety Frame for the medical and dental market. It is designed to be worn by itself to protect the user against specific laser radiation, or to be worn with Magnifying Loupes. Why do I need to wear laser eyewear?

The 'light' from medical and dental lasers can be concentrated and that light can cause permanent eye damage. It is possible to cause irreversible ocular injury with one direct or reflected laser beam. The laser light can be visible or invisible and the laser beam can be only as wide as a human hair.

- The frame design is light weight, wraps around the face, and has an removeable nose bridge. The intended market is Medical and Dental laser professionals.
- Lightweight, fit-over style
- Flexible earpiece for user comfort

**CE Scale Number:**

>190-315 D LB6 + I LB3 + R LB4 INVO CE  
 >315-1400 D LB6 + IR LB7 + M LB7Y INVO CE  
 >1400-11000 D LB3 + I LB4 INVO CE



The Renegade Safety Frame is identical to the PPE which is subject of the EU-type-examination certificate no C5877INVO/RO

DIN CERTCO Gesellschaft für  
 Konformitätsbewertung mbH  
 Alboinstr. 56  
 12103 BERLIN  
 GERMANY  
 Identification number of the Notified Body: 0196

The declaration of conformity is issued under the sole responsibility of the manufacturer

Signature:

Thomas Barrows  
 President

Date Issued: 10/9/2018