


PTCEIII – Blue Adjustable Goggle



- Laser and IPL protection
- CE Certified
- Includes silicone strap and hard case.
- Can be cold sterilized or autoclaved.
- Eye cups move independently.

Our patented PTCE Patient Goggle has been considered the top-of-the-line goggle for many years. Recent updates, to the nose bridge and eye cups, have provided the PTCE goggle with unmatched adjustability and ease-of-use. The nose bridge can swing 180 degrees and easily adjusts for different pupillary distances. The stainless-steel eyecups have been permanently coated with high grade silicone for added patient comfort.



Why is the PTCE SO Much Better?


Compare the old way to the Innovative Optics PTCEIII

How can you treat the inner canthus/
nose/glebella with this nose bridge?
You cannot move the brow bar.

VS

Our nose bridge EASILY moves out
of the way. Simply flip up and down
during treatment.

Competitors




They use gauze under their
goggle? Because their
goggle is uncomfortable.
This is a fire hazard you
should never use gauze
under a goggle.

VS

We have a Blue Silicone
Ring which you can
autoclave. This adds a layer
of comfort to the patient
and is not a fire hazard.

PTCE



They are using a **BLACK** Strap. You should
never use a black strap made of a material
you can't autoclave and disinfect.

VS

We have a **WHITE** silicone strap. It should
ALWAYS be white. You can autoclave the
strap or soak in disinfectant.

